

<b>Yipin Technology (Huizhou) Co., Ltd.</b>		Document Number	
		Version	A.0
<b>Technical Specifications</b>		Page	1/1
		File Name	HR100 Hair Removal Device Technical Specifications

1. Product Name: Hair Removal Device
2. Model: HR100
3. Performance Parameters;
  - 1) AC100-240V input, DC12V 4A output
  - 2) Irradiation Area: 3.0cm\*1.2cm = 3.6cm<sup>2</sup> or 2.0\*1.0cm = 2.0cm<sup>2</sup>
  - 3) Emission lifespan: 900,000 times
  - 4) Introduction to Touch Button Functions:
    - a) The middle touch button can be turned on/off by gently touching it with your finger for two seconds;
    - b) The left touch button allows you to adjust the gear level up or down by gently touching it when the device is turned on;
    - c) The right button can be toggled to the dual mode of light emission and cooling by gently touching it for 2 seconds. Another 2-second touch switches to the cold compress mode, and a further 2-second touch returns to the laser emission mode;
  - 5) Powering on normally:

Hold the machine handle, plug in the power adapter, observe if the top of the machine lights up with a green indicator. Then, gently touch the middle touch button to power on. Place the front aluminum surface in contact with the skin and observe if the indicator light is flashing rapidly. If it is flashing rapidly, press the button on the handle once to release a laser pulse. The device is ready to recharge when the indicator light flashes rapidly again, and you can press the handle button again to emit light;
  - 6) Automatic continuous light emission:

After pressing on the button on the handle for 3 seconds, release the button when the lightning icon in the bottom left corner of the screen lights up. At this point, as long as the front aluminum surface is in contact with the skin, it will automatically emit continuous flashes. Please wear protective goggles. In continuous emission mode, the green indicator light on the top of the machine will not flash rapidly
  - 7) Aging Mode:

Press and hold the button on the handle without releasing it, then plug in the power adapter. Immediately after, lightly touch the right touch button to activate the aging mode. In the aging mode, the first level emits light continuously 100 times, then jumps to the fifth level and emits light continuously 100 times. After the normal aging process is complete, the display screen will alternate between "---PS---" and the number "900,000". You can unplug the power directly;
  - 8) Product Dimensions: 198\*106\*73mm;
  - 9) Color Box Dimensions: 240\*210\*85mm;
  - 10) Carton Dimensions: 500\*440\*280mm (Each box contains 12 units) ;
  - 11) The net weight/gross weight for one unit is 0.337kg/0.944kg.  
The net weight/gross weight for each carton is 4.044kg/12.328kg;
4. Compliance with Design:
  - 1) ROHS Environmental Protection
  - 2) GB/T 2423.8-1995 Free Fall
  - 3) CE EMC EU Conducted Radiation
  - 4) FCC US Conducted Radiation



Designer: Song Qilin

Check:

Approval: